Directional Lead Positioner

Directions for Use

L011-90 (Rev C0, 2021-05-03)

Contains directions for the following products:

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English

Reference the microTargeting™ STar™ Drive System Directions for Use (L011-1007, included with STar™ Drive System shipment) for proper use instructions, including important safety precautions and warnings.


Indications for Use

The drive systems are intended to be used with commercially available stereotactic systems for neurosurgical procedures which require the accurate positioning of microelectrodes, stimulating electrodes, or other instruments in the brain or nervous system.

Intended Use

The microTargeting™ Directional Lead Positioner is intended for use by neurosurgeons in a standard operating room environment to support compatible leads at a specific distance during implantations.

Symbol Key

	WARNING/Caution, consult instructions for important cautionary information.		Medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, 98/79/EC and Medical Device Regulation (EU) 2017/745.
	Indicates the need for the user to consult the instructions for use.		Telephone Number
Rx Only	Caution - Federal law (USA) restricts this device to sale by or on the order of a physician.		Authorized Representative in the European Community.
	In reference to "RX Only" symbol; this applies to USA audiences only.		European Conformity. This device fully complies with Medical Device Regulation (EU) 2017/745 and legal responsibilities as a manufacturer are with FHC, Inc., 1201 Main Street, Bowdoin, ME 04287 USA.
	Indicates the batch code so that the batch or lot can be identified.		An item that is known to pose hazards in all MRI environments.
	Indicates the catalog number so that the medical device can be identified.		Indicates the temperature limits to which the medical device can be safely exposed.
	Indicates the serial number so that a specific medical device can be identified.		The date when the medical device was manufactured.
	Indicates the use of Radio-Frequency Identification (RFID).		
	Indicates Medical Device		

Handling and Storage

Storage: Store the Directional Depth Stop and Directional Lead Holder at normal room temperatures between -34°C (-29° F) and 57°C (135°F). Do not exceed 135°F for long-term storage.

Handling and use: Handle the Directional Depth Stop and Directional Lead Holder with extreme care. These components may be damaged if excessive force or incorrect handling occurs. Do not use any tool to assemble, use or disassemble the Directional Lead Positioner.

The products referenced in this document satisfy the FDA's Unique Device Identification Direct Marking requirement through the use of Radio-Frequency Identification (RFID). The RF air protocol is compliant with ISO-18000-6C EPC Global Class 1 Generation 2, and operates at frequencies in the band of 902-928MHz (UHF). The Unique Device Identification number retrieved from the RFID can be used to access product information through the FDA's Global Unique Device Identification Database (GUDID), which is available on the FDA's website.

Warnings and Cautions

WARNING: All parts are shipped non-sterile and must be sterilized before use. Follow sterilization protocols below.

WARNING: Components should be examined after each sterilization cycle for damage and function. Users should be aware that the effects of unvalidated sterilization protocols could result in damage to the components and affect their function or performance.

WARNING: When there is a lead in the brain, every effort should be made to minimize rotational forces to the tube or depth stop as it can translate into rotation movements of the lead in the brain.

Caution: Federal Law restricts this device to sale by or on the order of a physician.

Caution: The Directional Lead Positioner is not MRI safe.

Disposal: Dispose of accessories according to hospital protocol.

Cleaning Instructions

There are 2 methods to clean the Directional Lead Positioner. Select one.

Method: Manual		Container: Sterisuite Case	
Protocol			
Phase	Duration	Components/Notes	Detergent Type
Soak	5 minutes in detergent solution	Immerse all parts separated from each other.	Asepti Wash Plus Liquid
Wipe		Use detergent dampened cloth to wipe parts. Use brush to reach hard to clean areas	
Sonicate	10 minutes minimum	Tray fully loaded with parts in sonication unit with detergent	Asepti Wash Plus Liquid
Rinse		Reverse osmosis/de-ionized water	
Dry		Use clean soft cloth	
Or			
Method: Automated		Container: Sterisuite Case	
Protocol			
Phase	Recirculation Time	Water Temperature	Detergent Type
Pre-Wash 1	2 minutes	Cold tap water	N/A
Enzyme Wash	2 minutes	Hot tap water	Asepti Wash Plus or Sekusept AR
Wash 1	2 minutes	65.5°C	Asepti Wash Plus or Sekusept AR
Rinse 1	2 minutes	Heated tap water	N/A
Thermal Disinfection (optional)	5 minutes	90°C	N/A
Pure Water Rinse	10 seconds	Heated	Asepti Rinse or Sekusept FNZ or Sekumatic Multiclean
Dry	7 minutes	115°C	N/A

Listed are the detergents and cycles that FHC has validated. Detergents listed are from Ecolab. If other neutral or alkaline detergents are used, testing should be done by the hospital to ensure product is not damaged. Detergents should be prepared per manufacturer's recommendations.

Sterilization Instructions

Load items to be sterilized into the Sterisuite Case. For proper placement refer to STar™ Drive system DFU L011-1007. Sterilize using one of these validated protocols:

Prevacuum (Steam):	
Temperature:	132°C (270°F)
Preconditioning Pulses:	3
Exposure Time:	4 minutes
Minimum Dry Time:	60 minutes
Cycle was validated using Halyard Health H400 wrap	
Temperature:	134°C (273°F)
Preconditioning Pulses:	3
Exposure Time:	18 minutes
Minimum Dry Time:	60 minutes
Cycle was validated using Halyard Health H300 wrap	
Gravity (Steam):	
Temperature:	132°C (270°F)
Exposure Time:	30 minutes
Minimum Dry Time:	60 minutes
Cycle was validated using Halyard Health H300 wrap	
Double Pouched in self-sealing pouches	
Prevacuum (Steam):	
Temperature:	134°C (273°F)
Preconditioning Pulses:	4
Exposure Time:	18 minutes
Minimum Dry Time:	20 minutes
Cycle was validated using Cardinal Health self-sealing pouch CAT #92713 and #92152	

Maintenance and Repair

Scheduled maintenance: The Directional Depth Stop and Directional Lead Holder should be inspected visually prior to each use for physical damage, or poor fit due to wear or residue buildup. If any component shows wear or damage that could interfere with proper function, please contact FHC for repair or replacement. None of the components require lubrication of any kind.

Repair: All FHC products are unconditionally guaranteed against defects in workmanship for one year from the date of shipment provided they have been exposed to normal and proper use. Should service or repair be required, please contact FHC for return instructions at 1-800-326-2905 (US & Canada) or +1-207-666-8190.

Pre-Use Check


The following steps should be taken prior to performing the actual procedure:

1. Check the integrity of the plastic tip on the depth stop screw.
2. Confirm the depth stop cannot be rotated in the lead holder when fully inserted
3. Test mount the lead holder on the STar™ Drive carriage and confirm there is no looseness or any stripped threads.
4. Inspect Directional Lead Holder and confirm pin is present and not damaged.

Specifications

Dimensions

- Discrete Rotational Increments: 30°
- Maximum of 120° of rotation in either direction from the 0° point
- Distance to target dimensions are equivalent to existing products (L011-1007)
- Compatible with Abbot and Boston Scientific Directional leads.


Material Profile

Depth Stop Screw:	316 Stainless Steel w/ PEEK tip
Depth Stop Body:	Natural Radel
Lead Holder:	Aluminum 6061-Hard coat anodize black with PTFE
Lead Holder Thumbknobs:	316 Stainless Steel
Lead Holder Pin:	316 Stainless Steel


Procedure

DBS Lead Preparation (should be done in a sterile area)


1. Attach Directional Lead Holder to STar™ Drive Lead Measurement Fixture then insert the Directional Depth Stop into holder, positioned at zero, then tighten lead holder depth stop locking thumbknob.


- Loosen depth stop screw and insert lead through the Directional Depth Stop until tip of the lead is at the proper position on the measurement fixture. Orient the lead so that the alignment indicator is facing up and aligned with the depth stop screw. Lightly tighten the depth stop screw to secure the lead.

- Loosen the lead holder depth stop locking thumbknob then lift and rotate the Directional Depth Stop to the desired angle depending on the Star™ drive orientation on the stereotactic system. Re-insert the Directional Depth Stop and secure using the thumbknob.

4. Remove Directional Lead Holder (with Directional Depth Stop and lead attached) from measurement fixture.

DBS lead Insertion

1. The lead insertion and holder attachment procedure is identical to the standard lead holder and depth stop (refer to L011-1007).

Securing the DBS Lead

1. The cannula retraction and lead securing procedure is identical to the standard lead holder and depth stop (refer to L011-1007)

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