microTargeting™ DP Stylets

Directions For Use

L011-81-01 (Rev B0, 2021-05-24)

Contains directions for the following products:

66-ST-x (where x represents the following variations- Length: 15mm – 300mm, Outside diameter: 0.5-2.0mm, Collar size: 1-5mm diameter, "P" indicates presterile)

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Indications for use

The FHC microTargeting™ DP Stylets are nonpowered hand-held instruments intended to cut, hold or manipulate tissue during stereotactic functional neurosurgical procedures.

Intended use

The DP Stylets are a part of the microTargeting SStar Drive System and is intended for use by neurosurgeons in a standard operating room environment to position and move devices relative to a planned stereotactic system trajectory and target.

Symbol Key

	WARNING /Caution, consult instructions for important cautionary information.		Medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, 98/79/EC and Medical Device Regulation (EU) 2017/745.
	Consult instructions for use.		Telephone number
	In reference to "Rx only" symbol; this applies to USA audiences only.		Authorized Representative in the European Community.
Rx Only	Caution- Federal law (USA) restricts this device to sale by or on the order of a physician.		European Conformity. This device fully complies with MDD Directive 93/42/EEC and legal responsibilities as a manufacturer are with FHC, Inc., 1201 Main Street, Bowdoin, ME 04287 USA.
	Indicates a Medical Device		0413
	Do not re-use; intended for one use on a single patient, during a single procedure.		Do not use if package is damaged or opened.
	Medical device that is not to be re-sterilized.		Date after which the medical device is not to be used.
	Non-pyrogenic		Indicates the temperature limits to which the medical device can be safely exposed.
	Indicates the catalog number so that the medical device can be identified.		Indicates the range of humidity to which the medical device can be safely exposed.
	Indicates the batch code so that the batch or lot can be identified.		Sterilized using ethylene oxide.
	The date when the medical device was manufactured.		Medical device that has not been subjected to sterilization process.
	Instructions for end of life disposal.		

microTargeting™ Is a trademark of FHC, Inc.

Sterile DP Stylets

STERILE EO

WARNINGS

- Sterile Medical Device – Do NOT re-sterilize.
- Do not use the contents if there is any evidence of damage to the package or package seal that could compromise sterility.
- For single patient use only
- Do not reuse; reusing single-use medical devices could lead to serious patient injury.

Non Sterile DP Stylets

NON STERILE

- If your microTargeting DP Stylets were shipped in a standard plastic box, DO NOT attempt to sterilize the stylets using this package. Remove microTargeting DP Stylets from the package prior to sterilization. FHC recommends use of the FHC sterilization tray (holds seven components).
- FHC has validated and recommends the following steam sterilization parameters:

Prevacuum wrapped (in 2 layers of 1-ply polypropylene wrap ⁽¹⁾) Preconditioning pulses: 3 Exposure time: 4 minutes at 132°C (270°F) Minimum dry time: 30 minutes ⁽¹⁾ Cycle was validated using Halyard Health H600 wrap	Gravity unwrapped Exposure time: 10 minutes at 132°C (270°F)
--	--
- FHC has validated and recommends the following manual cleaning protocol:
 - Immediately following use, thoroughly rinse each stylet under tap water.
 - Immerse all components in enzymatic wash solution for 1 minute. Remove from the wash solution, rinse with tap water, then distilled water.
 - Dry devices using a clean soft cloth.
- There is no restriction on the number of reuses providing user performs pre-use check.

Safety Information

- Not intended for implantation

Contraindications

- microTargeting DP Stylets are not suited for chronic implantation.

WARNINGS

Rx Only • **Rx only: Caution-** Federal law (USA) restricts this device to sale by or on the order of a physician.

- In the event that a patient must be defibrillated via electric shock, the microTargeting DP Stylet should be withdrawn prior to defibrillation.
- To prevent electrical shock hazard DO NOT connect microTargeting DP Stylets to any line voltage source or any unknown power source.
- When there is a stylet in the brain, every effort should be made to minimize lateral forces to the drive/stereotactic system.

Pre-Use Check

WARNING: Delicate stylets should always be carefully inspected before use because damage to them can impact patient safety. Handle stylets carefully to prevent bending.


WARNING: Check to be certain that stylets are straight; with the collar over the edge, roll the microTargeting™ DP Stylet over a flat surface or in the measuring fixture slot.

Directions For Use

1. Adjust the stereotactic positioning system for the desired trajectory and depth.
2. Insert the stylet into the stereotactic positioning system.
3. Puncture the dura.
4. After use, withdraw the stylet from the stereotactic positioning system.
5. Refer to the directions for use associated with the insertion tubes and microelectrodes.

6. Proceed with cleaning instructions for non sterile microTargeting DP Stylets, otherwise place the stylet in an approved sharps disposal container according to hospital protocol.

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