



microTargeting™ Solution for PMT® Depth Electrode Placement

Patient-specific stereotactic system that allows a surgical team to quickly place an unlimited number* of PMT® SEEG depth electrodes using PMT® skull bolts, drill and wrench

FHC's microTargeting solution for PMT® SEEG depth electrode placement provides guidance for the PMT® 2.4mm drill and Anchor Bolt Driver to accurately place PMT® Anchor Bolts and Depth Electrodes using WayPoint™ bone anchors serving as fiducials, and a patient-specific microTargeting Multi-Oblique Platform created with WayPoint Navigator software. Over twenty independent trajectories may be configured with a single platform. Once the platform is attached, all of the trajectories are available for implantation. This allows for rapid progression between trajectories as well as the possibility of simultaneous use by multiple surgeons.


FHC also offers both a microTargeting solution for Visualase® or Monteris® Laser Ablation and one for NeuroPace RNS® System, that allow both localization and treatment to be performed using the same technology.

* "There is no practical limitation to the number of electrodes that can be implanted."

Dewan MC et al. J Neurosurg. 2017 Dec 15:1-9. doi: 10.3171/2017.6.JNS17881

microTargeting™ and WayPoint™ are trademarks of FHC, Inc.

- ✓ WayPoint bone anchors serving as fiducials are placed in clinical setting 7-28 days before surgery
- ✓ Patient-specific platform quickly attaches to bone anchors in OR
- ✓ No restrictions on positioning of the anchors, selecting the trajectory entry point or orientation
- ✓ Anchors and scans may be reused for subsequent microTargeting treatment solutions
- ✓ Supports implantation of all SEEG depth electrodes using a single platform
- ✓ Platform delivered to hospital in 5 days



Procedure Overview:

1. Use 2.4mm Tool Guide (66-MO-04) to guide PMT® drill through Skull.
2. Use 2.5mm Tool Guide (66-MO-05) to screw PMT® bolt into the hole using PMT® wrench.
3. Use PMT® drill, PMT® drill stop and ruler to determine distance from top of bolt to target (refer to figure) and record on trajectory depth sheet.
4. Remove guide and place solid PMT® cap on Bolt to prevent CSF leakage.
5. Repeat steps 1-5 for all trajectories, then remove platform.
6. For each trajectory, position bolt cap on lead to distance in step 3 then insert and tighten cap.

Ordering Information:

Reusable tools/planning software:

66-ST-MO – Sterilization tray, ruler, cleaning brush
 66-MO-04 – 2.4mm tool Guide
 66-MO-05 – 2.5mm tool Guide
 66-WP-NV (4.5) – WayPoint Navigator

Per case:

MP-KIT-P-EO – Multi trajectory oblique platform

Additional Items Required for Operation:

(vendor) PMT® 2.4mm Drill
 (vendor) PMT® Anchor Bolt
 (vendor) PMT® Anchor Bolt Driver
 (vendor) PMT® depth electrodes

L022-08 (Rev C2, 2020-06-15)



FHC, Inc.
 1201 Main Street
 Bowdoin, ME 04287 USA
 Fax: +1-207-666-8292
 www.fh-co.com



FHC Europe
 (TERMOBIT PROD srl)
 42A Barbu Vacarescu Str, 3rd Fl
 Bucharest 020281 Sector 2
 Romania

FHC Latin America
 Calle 6 Sur Cra 43 A-200
 Edificio LUGO Oficina 1406
 Medellín-Colombia