Distributed by

Medtronic

UNIBODY™

Bone Fiducials

Technical Manual

L011-500-03 (Rev A0, 2019-08-27)

Contains Directions for the following products:
FM-4007, FM-4010, FM-4013

FHC, Inc.
1201 Main Street
Bowdoin, ME 04287 USA
www.fh-co.com

Distributed By:
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604

EC REP

FHC Europe
(TERMOBIT PROD srl)
42A Barbu Vacarescu Str, 3rd Fl
Bucharest 020281 Sector 2
Romania

Table of contents

English	4
Français	8
Deutsch	12
Italiano	17
Español	21
Português	25
Nederlands	29
Svenska	33
Dansk	37
Norsk	41
Čeština	45
Magyar	49
Slovenčina	53
Ελληνικά	57
Türkçe	61
Polski	65
Русский	69

English

Explanation of symbols on package labeling

WARNING/Caution, consult instructions for important cautionary information.

Consult instructions for use.

European Conformity. This device fully complies with MDD Directive 93/42/EEC and legal responsibilities as a manufacturer are with FHC, Inc., 1201 Main Street, Bowdoin, ME 04287 USA.

Authorized Representative in the European Community.

Medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC.

Rx Only

Caution- Federal law (USA) restricts this device to sale by or on the order of a physician.

In reference to "Rx only" symbol; this applies to USA audiences only.

Indicates the batch code so that the batch or lot can be identified.

Indicates the catalog number so that the medical device can be identified.

Indicates a medical device that has been sterilized using ethylene oxide.

Indicates the date after which the medical device is not to be used.

Do not re-use; intended for one use on a single patient, during a single procedure.

Indicates a medical device that is not to be resterilized.

Do not use if package is damaged or opened.

Indicates the temperature limits to which the medical device can be safely exposed.

Medtronic® is a registered trademark of Medtronic, Inc.

Unibody™ is a trademark of Medtronic, Inc.

Device description

Unibody Models FM-4007, FM-4010, and FM-4013 Bone Fiducials provide fixed reference points during neurosurgical procedures.

Package contents

- Bone fiducials:
 - 7-mm bone fiducials (6) (FM-4007)
 - 10-mm bone fiducials (6) (FM-4010)
 - 13-mm bone fiducials (6) (FM-4013)
- Protective caps (6)
- Screwdriver guides (6)
- Product literature

Device specifications

Table 1. Material of components in the Models FM-4007, FM-4010, and FM-4013 packages

Component	Material	Material contacts human tissue
Bone fiducial	Titanium	Yes
Protective cap	Thermoplastic	Yes
Screwdriver guide	Polycarbonate	No

Indications for Use

The Unibody fiducials are part of the Nexframe Stereotactic System which is intended to provide stereotactic guidance for the placement and operation of instruments or devices during planning and operation of neurological procedures performed in conjunction with the use of an image Guided Workstation System using preoperative MR and/or CT imaging. These procedures include biopsies, catheter placement and electrode introduction.

Intended Use

The Unibody FM-4007, FM-4010, and FM-4013 Bone Fiducials are intended for use by neurosurgeons as fixed reference points for presurgery CT imaging for patients requiring stereotactic surgery in standard operating room environment.

WARNINGS

 **Single use only** – This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Component packaging – Do not use a component if the following circumstances have occurred:

- The storage package has been pierced or altered because component sterility cannot be guaranteed and infection may occur.
- The component shows signs of damage because the component may not function properly.
- The use-by date has expired because component sterility cannot be guaranteed and infection may occur.

Instructions for use

WARNINGS:

- The bone fiducials are not intended for permanent implant.
- Do not place the bone fiducials in a cranial suture or in cranial bone that is damaged or diseased, or is less than 5 mm thick. If the bone fiducials are placed in diseased or unstable bone, the following can result:
 - The bone fiducials will not be properly anchored, causing potential fiducial movement and registration errors.
 - The bone fiducials will be unstable, increasing the risk of infection, damage to the dura, possible CSF leakage, and damage to the neural (brain) tissue.

Placing the bone fiducials

1. Prep and incise the patient at the desired location for each bone fiducial placement. A minimum of four locations should be chosen. The locations should cover the entire cranium and be noncoplanar.
2. Insert a bone fiducial through a screwdriver guide sleeve.

3. Position and hold the guide sleeve to ensure the bone fiducial is placed orthogonally to the skull surface.
4. Using a powered screwdriver with a 1.6-mm cross-type blade, place the blade through the guide sleeve hole and drive the bone fiducial until it is flush with the skull.
5. Repeat steps 2-4 for the remaining bone fiducial placements.
6. Use a manual screwdriver with a 1.6- mm cross-type blade to ensure each bone fiducial is fully seated.

Caution: To avoid stripping, do not overtighten the bone fiducials.

7. Insert a protective cap over each bone fiducial.
8. Perform CT imaging and acquire the patient's CT scans.

Caution: Keep the patient's head immobile during imaging.

Preoperative use

1. Remove the protective caps and drape the patient.
2. According to the image-guided surgery (IGS) system instructions:
 - a. Select the center of each bone fiducial.
 - b. Begin the registration process.
 - c. Touch the center of each bone fiducial divot.

Caution: Be sure to touch the bottom of the conical center of each bone fiducial divot to locate the physical center of each fiducial location.

3. Complete the surgery.

Postoperative removal

1. Remove all bone fiducials and discard in accordance with hospital protocol.
2. Clean and close all fiducial incisions.

Disposal Instructions

Dispose of the fiducials according to Hospital Protocol.

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