

Guideline 4000 LP+™

Neuromodulation Targeting System

microTargeting™

Cart and case options for set-up

Portable and highly configurable, the LP+ components can be arranged in any number of ways to meet the specific needs of the surgical team.



Ordering Information

MT-LPP-MEPK

Mobile electrophysiology kit - includes Guideline 4000 LP+, mT Controller, single case, QuickWave, serial-USB cable, and hard case

MT-LPP-CART

Cart option, price includes installation (US only)

MT-LPP-HCSE

Single shipping case, holds LP+, laptop, mT Controller, remotes, and speaker (airplane check)

MT-LPP-SCSE

Two shipping cases, one for the LP+, one for other items (airplane carry-on)

L012-LP+ Options RevB

Cart option – The entire LP+ system can be arranged on a wheeled cart for permanent installation. Includes speaker mounts, mouse pad, and a locking cabinet for equipment storage.

Hard case option – All components fit in one foam-padded, lockable, wheeled transport case with telescoping handle. Suitable for airplane check, car transport, or shipment.



Soft case option – The LP+ and laptop fit in a small, soft case. A second small case carries accessories and associated equipment. Suitable for airplane carry-on.

FHC



FHC, Inc.
1201 Main Street
Bowdoin, ME 04287
Fax: (207) 666-8292
Email: fhcinc@fh-co.com
www.fh-co.com



24 hour technical service:
1-800-326-2905 (US & Can)
(207) 666-8190

EC REP

FHC Europe
(TERMOBIT PROD srl)
129 Barbu Vacarescu Str, Sector 2
Bucharest 020272
Romania

"Innovation through collaboration"