

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.;
Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Electrodes

The following section has been included to assist you in calculating the price of the various types of microelectrodes we carry. Please refer to the catalog pages for unique specifications available for each type.

UE1234567890n/n/n Metal Microelectrodes

FHC's UE microelectrodes are sold in packages of 12 with the exception of Platinum/Iridium, which is sold in packages of 3. Prices are per electrode and subject to change without notice. To calculate the cost of metal microelectrodes, follow the formula noted, substituting the price shown for each specification.

Formula: $1 \times 2 \times 3 + (4+5+6+7+8+9+10) = \text{Price per microelectrode}$

1 Material

W & S: \$15.00
P: \$28.00

2 Electrode Length

Multiply the Material Type by the factor listed below:

S & M: x 1.0
L: x 1.25
V: x 3.0

3 Shank Diameter

Take the sum from above and multiply by the factor listed below:

C: x 1.50
D: x 1.30
E: x 1.15
F & G: x 1.00
H: x 1.10

4 Final Taper Angle

A,C,E,G & J: No charge
B,D,F,H & K: Add \$0.55ea.

5 Microelectrode Profile

S,V & C: No extra charge
T & L: Add \$1.10ea.

6 Insulation Options

E & N: No charge
G & S: Add \$8.50
M: Add \$1.10

7 Impedance

Add \$0.50ea. if ordering more than one impedance per package

8 Tip Conditioning

N: No charge
P & T: Add \$2.00

Example of a standard length tungsten microelectrode:

Catalog Number: **UEWMFGSEDNNM**
Calculation: $(\$15.00) \times (1.0) \times (1.0) + (0+0+0+0+0+0+\$1.10)$
Price per electrode: \$16.10 or \$193.20 per package of 12

Example of a short length platinum/iridium microelectrode:

Catalog Number: **UEPSECSGDN5H**
Calculation: $(\$28.00) \times (1.0) \times (1.15) + (0+0+8.50+0+0+4.75+0)$
Price per electrode: \$45.45 or \$136.35 per package of 3

9 Shank Modifications

N & 1: No charge
2: Add \$0.60
3 up to 150mm: Add \$5.50
Over 150mm: Add \$12.50
4 up to 150mm: Add \$3.40
Over 150mm: Add \$10.50
5: Add \$4.75
6: Add \$11.50

10 Termination Options

M & F: Add \$1.10
C & L: Add \$3.30
E,G,H & J: No charge

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Electrodes

CE12345n/n/n Cluster Electrodes

1 Number of Contacts

- 2: 2 Conductors \$66.00ea.
- 3: 3 Conductors \$70.00ea.
- 4: 4 Conductors \$75.00ea.
- 5: 5 Conductors \$80.00ea.

2 Electrode Material

- Multiply by the factor shown:
- B: x 2.0
 - C: x 1.6
 - D: x 1.4
 - E: x 1.6
 - F: x 1.4

Cluster electrodes are sold in packages of 3

To calculate the cost, follow the formula noted, substituting the price shown for each specification.

Formula: 1 x 2 = Price per Cluster Electrode

Example of a 2 Conductor Cluster Electrode:

Catalog Number: **CE2C75**
Calculation: (\$66.00) x (1.6)
Price per electrode: \$105.60 or \$316.80 per package of 3
(Please add \$0.10 for every mm over 100mm)

CB123456n/n/n Concentric Bipolar Electrodes

1 Outer Pole Diameter

- A: \$130.00ea.
- B: \$66.00ea.
- C: \$62.00ea.
- D: \$62.00ea.
- E: \$62.00ea.
- F: \$62.00ea.
- G: \$59.00ea.
- H: \$59.00ea.
- J: \$59.00ea.

3 Inner Pole Diameter

- B, C, E, G, J & L:
Multiply the price by 1.25
- D, F, H, K & M:
(Stainless Steel)
Multiply the price by 1.15

"M" Tip Configuration:
Add price of UE Metal Microelectrode

4,5,6 Length

(In millimeters) add \$0.10 for every millimeter over 100mm

Concentric Bipolars are sold in packages of 3

To calculate the cost, follow the formula noted, substituting the price shown for each specification.

Formula: (1 x 3) + 4,5,6 = Price per Concentric Bipolar Electrode

Example of a 125µm Outer Pole/25µm Pt/Ir Inner Pole CB Electrode:

Catalog Number: **CBARC75**
Calculation: (\$130.00) x (1.25) + (0)
Price per electrode: \$162.50 or \$487.50 per package of 3

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Electrodes

UN123456(ii) Microneurography Needles

1 Configuration

A: \$180.00/pkg 12
R: \$88.00/pkg 12
P: \$135.00/pkg 6 pair

2,3 Length

Up to 60mm: No Charge
60 to 130mm: Add \$12.00/pkg

4 Diameter

E: Add \$12.00/pkg of 12
F & G: No Charge

5 Impedance

No Charge

6 Termination


D,E,G,H & J: No Charge
R & S: \$30.00/pkg 12
T & L: \$44.00/pkg 12
U: \$65.00/pkg 12
M: \$15.00/pkg 12

Microneurography needles are packaged in boxes of 12 Active Needles, 12 Reference Needles, or pairs of 6 Reference and 6 Active Needles.

To calculate the price, follow the formula noted, substituting the price shown for each specification.

Formula: 1 + 2,3 + 4 + 6 = Package price of Microneurography Needles

Example of a Microneurography Active Needle:



Catalog Number: **UNA35F2U**
Calculation: (\$180.00) + (0) + (0) + (\$65.00)
Package Price for 12 Active Electrodes: \$245.00

MX12345(ii) Matrix Electrode

1,2 Number of Contacts + 1

3 Spacing

Multiply by 1.0 for 240 micron and greater spacing
Multiply by 1.2 for 220 micron and smaller spacing

4 Connectors

For Connector C or D add \$3.00 per contact


5 Material Cost of Individual Electrodes Used

Tungsten or Stainless Steel: \$18.50
Platinum/Iridium: \$45.00

For contacts greater than 9, please call for pricing. To calculate the cost of Matrix Electrodes, follow the formula noted, substituting the price shown for each specification.

Formula: (1 x 2 + 1) x (3) x (5) + (4) = Price per Matrix Electrode

Example of a Matrix Electrode:



Catalog Number: **MX21AES(DH1)**
Calculation: (1 x 2 + 1) x (1.2) x (18.50)
Price per MX electrode: \$66.60 or \$199.80 per package of 3

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.;
Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price
List L017W-02 US
Rev. 11/28/11

Electrodes

PB1234567n/n/n Parallel Bipolar Electrode

1 Base Price Configuration

All: \$60.00 each

2 Material


A: \$0.065/mm

B: \$0.11/mm

C & D: \$0.035/mm

Material Costs:

Total material charge is determined by multiplying the material cost by the length in mm as specified by factors 5, 6 and 7 on the tech sheet.



Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.;
Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Micropipettes

Series 10, 11, 12, 13, 14, or 17-SD-TL

| | | | | | |
|----------------|-----------|---------------------------------|----------------|-----------|-------------|
| Specify | SD | 10 = 1.0mm | Specify | TL | S = 5-7mm |
| (Shank Diam.) | | 12 = 1.2mm | (Taper Length) | | M = 10-12mm |
| | | 15 = 1.5mm | | | L = 15-17mm |
| | | 20 = 2.0mm (single barrel only) | | | |

e.g. 12-12-S = 2 barrel; ea barrel shank 1.2mm OD, 7mm taper length

*An additional charge of \$2.50 each (1-3 barrel) or \$2.75 each (4-7 barrel) will be added for specific tip modification requirements. Please contact our Technical Department for details.

Prices listed are for each piece; the package quantity is listed in the description.

| Item | Description/Package Quantity | Price Ea. |
|----------|--|-----------|
| 10-SD-TL | 1 Bar .12% Wall/Borosil. Micropipette (pk.48) | \$ 6.42 |
| 11-SD-TL | 1 Bar .125% Wall/Borosil. Micropipette (pk.48) | \$ 6.42 |
| 12-SD-TL | 2 Bar .125% Wall/Borosil. Micropipette (pk.48) | \$ 7.06 |
| 13-SD-TL | 3 Bar .125% Wall/Borosil. Micropipette (pk.48) | \$ 7.57 |
| 14-SD-TL | 4 Bar .125% Wall/Borosil. Micropipette (pk.48) | \$ 7.91 |
| 17-SD-TL | 7 Bar .125% Wall/Borosil. Micropipette (pk.48) | \$ 9.05 |
| 18-20-1 | Ag/AgCl Electrode | \$ 48.40 |
| 18-20-2 | Ag/Ag/Cl Electrode/Platinum Blacked | \$ 52.80 |
| 19-20-1 | Micropipette Filling Kit | \$ 99.00 |
| 19-20-2 | 30 ga x 3" Filling Needle (pk.5) | \$ 17.60 |
| 19-20-3 | Micropipette Filling Jar | \$ 44.00 |

Capillary Tubing

Prices listed are for each piece; the package quantity is listed in the description. Standard length is 100mm (4"). Many types also stocked in 75mm (3") lengths. Other lengths also available; package quantity is adjusted for different lengths.

| Item | Description/Package Quantity | Price Ea. |
|----------|--|-----------|
| 27-30-0 | Borosil 1.0mm OD x 0.75mm ID 100mm, 250/Pack | \$ 0.35 |
| 27-30-1 | Borosil 1.0mm OD x 0.5mm ID 100mm, 250/Pack | \$ 0.35 |
| 27-31-0 | Borosil 1.2mm OD x 0.9mm ID 100mm, 250/Pack | \$ 0.35 |
| 27-31-1 | Borosil 1.2mm OD x 0.6mm ID 100mm, 250/Pack | \$ 0.35 |
| 27-32-0 | Borosil 1.5mm OD x 1.12mm ID 100mm, 250/Pack | \$ 0.55 |
| 27-32-1 | Borosil 1.5mm OD x 0.75mm ID 100mm, 250/Pack | \$ 0.55 |
| 27-33-1 | Borosil 2.0mm OD x 1.0mm ID 100mm, 250/Pack | \$ 0.59 |
| 30-30-0 | Borosil 1.0mm x 0.75mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-30-0X | Borosil 1.0mm OD x 0.85mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-30-1 | Borosil 1.0mm x 0.5mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-30-2 | 2 Barrel Borosil 1.0mm OD with Fiber, 100mm, 100/Pack | \$ 0.57 |
| 30-30-3 | 3 Barrel Borosil 1.0mm OD with Fiber, 100mm, 100/Pack | \$ 0.70 |
| 30-30-4 | 4 Barrel Borosil 1.0mm OD with Fiber, 100mm, 100/Pack | \$ 0.85 |
| 30-30-7 | 7 Barrel Borosil 1.0mm OD with Fiber, 100mm, 100/Pack | \$ 1.25 |
| 30-31-0 | Borosil 1.2mm x 0.9mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-31-1 | Borosil 1.2mm OD x 0.6mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Capillary Tubing continued

Prices listed are for each piece; the package quantity is listed in the description.

| Item | Description/Package Quantity | Price Ea. |
|----------|--|-----------|
| 30-31-1T | Borosil 1.2mm OD x 0.4mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-31-2 | 2 Barrel Borosil 1.2mm OD with Fiber, 100mm, 100/Pack | \$ 0.57 |
| 30-31-3 | 3 Barrel Borosil 1.2mm OD with Fiber, 100mm, 100/Pack | \$ 0.70 |
| 30-31-4 | 4 Barrel Borosil 1.2mm OD with Fiber, 100mm, 100/Pack | \$ 0.84 |
| 30-31-7 | 7 Barrel Borosil 1.2mm OD with Fiber, 100mm, 100/Pack | \$ 1.25 |
| 30-32-0 | Borosil 1.5mm OD x 1.12mm ID with Fiber, 100mm, 250/Pack | \$ 0.55 |
| 30-32-1 | Borosil 1.5mm OD x 0.75mm ID with Fiber, 100mm, 250/Pack | \$ 0.55 |
| 30-32-2 | 2 Barrel Borosil 1.5mm OD with Fiber, 100mm, 100/Pack | \$ 0.85 |
| 30-32-3 | 3 Barrel Borosil 1.5mm OD with Fiber, 100mm, 100/Pack | \$ 1.00 |
| 30-32-4 | 4 Barrel Borosil 1.5mm OD with Fiber, 100mm, 100/Pack | \$ 1.08 |
| 30-32-7 | 7 Barrel Borosil 1.5mm OD with Fiber, 100mm, 100/Pack | \$ 1.33 |
| 30-33-0 | Borosil 2.0mm OD x 1.5mm ID with Fiber, 100mm, 250/Pack | \$ 0.59 |
| 30-33-1 | Borosil 2.0mm OD x 1.0mm ID with Fiber, 100mm, 250/Pack | \$ 0.59 |
| 30-33-1T | Borosil 2.0mm OD x 0.7mm ID with Fiber, 100mm, 250/Pack | \$ 0.59 |
| 30-34-2 | 2 Barrel Iontophoresis Tubing 100mm, 100/Pack | \$ 0.84 |
| 30-50-2 | Borosil 0.5mm OD x 0.25mm ID 100mm, 250/Pack | \$ 0.35 |
| 30-50-3 | Borosil 0.5mm OD x .25mm ID with Fiber, 100mm, 250/Pack | \$ 0.35 |
| 30-50-6 | Borosil 0.7mm OD x 0.35mm ID/ | \$ 0.35 |
| 30-50-8 | Borosil 0.7mm OD x 0.35 mm ID 100mm, 250/Pack | \$ 0.35 |

Electrode Fabrication Supplies

Prices listed are for each piece; the package quantity is listed in the description.

| Item | Description/Package Quantity | Price Ea. |
|-------------|---|-----------|
| 30-35-1 | Stainless Steel Rods .005"x 6" 250/pack | \$ 0.57 |
| 30-35-2 | Stainless Steel Rods .010"x 6" 250/Pack | \$ 0.61 |
| 30-35-3 | Stainless Steel Rods .020"x 6" 250/Pack | \$ 0.64 |
| 30-36-1-075 | Tungsten Rods .005" x 3" 250/Pack | \$ 1.32 |
| 30-36-1B | Tungsten Rods .005" x 7" 250/Pack | \$ 2.21 |
| 30-36-1B-02 | Tungsten Rods .003" X 7" 250/Pack | \$ 2.21 |
| 30-36-2-075 | Tungsten Rods .010" x 3" 250/Pack | \$ 1.32 |
| 30-36-2B | Tungsten Rods .010" x 7" 250/Pack | \$ 2.21 |
| 30-36-8-075 | Tungsten Rods .008" x 3" 250/Pack | \$ 1.32 |
| 30-36-8B | Tungsten Rods .008" x 7" 250/Pack | \$ 2.21 |
| 30-37-1 | Polyimide Tubing .0095"x.008" 10 feet / Pack | \$ 17.85 |
| 30-37-1-01 | Polyimide Tubing .0086"x.0071" 10 feet / Pack | \$ 17.85 |
| 30-37-1-02 | Polyimide Tubing .0083"x.0063" 10 feet / Pack | \$ 17.85 |
| 30-37-2 | Polyimide Tubing.0133"x .0113" 10 feet / Pack | \$ 17.85 |
| 30-37-3 | Polyimide Tubing.0162"x.0142" 10 feet / Pack | \$ 17.85 |
| 30-37-6 | Polyimide Tubing.0199"x.0179" 10 feet / Pack | \$ 17.85 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.;
Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Electrode Fabrication Supplies continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|---------|---|-----------|
| 30-37-7 | Polyimide Tubing .0246"x.0226" 10 feet / Pack | \$ 17.85 |
| 30-37-8 | Polyimide Tubing .0304"x.0284" 10 feet / Pack | \$ 17.85 |
| 30-38-1 | Pt/Ir Rods .005" x 2.5" 50 / Pack | \$ 7.01 |
| 30-38-2 | Pt/Ir Rods .010" x 2.5" 50 / Pack | \$ 14.20 |
| 30-39-1 | SS Fine Wire .001" (25micron) Insulated, 10 feet / Pack | \$ 5.63 |
| 30-39-2 | SS Fine Wire .001" (25micron) Uninsulated, 10 feet / Pack | \$ 5.63 |
| 30-39-3 | SS Fine Wire .002" (50micron) Insulated, 10 feet / Pack | \$ 5.63 |
| 30-39-4 | SS Fine Wire .002" (50micron) Uninsulated, 10 feet / Pack | \$ 5.63 |
| 30-40-1 | Female Amphenol Pins 100/Pack | \$ 1.58 |
| 30-40-2 | Male Amphenol Pins 100/Pack | \$ 1.60 |
| 30-40-4 | 53-pos. strip for male pins | \$ 25.20 |
| 30-40-5 | 53-pos. strip for female pins | \$ 25.20 |
| 30-40-8 | Crimping Tool | \$ 493.50 |
| 30-40-9 | Die for Amphenol Pins | \$ 115.50 |

Devices and Accessories

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|----------|--|-----------|
| 40-10-1 | Rack Mount Assembly | \$ 105.00 |
| 40-10-2 | 3-Module Cabinet | \$ 237.00 |
| 40-10-3 | Blank Front Panel (A) 7x17.8cm | \$ 28.00 |
| 40-10-4 | 6 Module Cabinet | \$ 325.00 |
| 40-10-5 | Blank Module (A) 7x17.8x22.8cm | \$ 94.00 |
| 40-10-6 | Blank Module (C) 14x17.8x22.8cm | \$ 121.00 |
| 40-10-7 | Blank Module (D) 21x17.8x22.8cm | \$ 144.00 |
| 40-10-8 | Blank Front Panel (C) 14x17.8cm | \$ 32.00 |
| 40-10-9 | Blank Front Panel (D) 21x17.8cm | \$ 34.00 |
| 40-11-1 | Banana/Pin Tip Adaptor - Black 5/Pack | \$ 7.50 |
| 40-11-10 | Alligator Clip smooth | \$ 3.25 |
| 40-11-2 | Banana/Pin Tip Adaptor - Red 5/Pack | \$ 6.50 |
| 40-11-3 | Binding Post/BNC Adaptor Black 5/Pack | \$ 5.25 |
| 40-11-4 | Binding Post/BNC Adaptor Red 5/Pack | \$ 5.25 |
| 40-11-5 | BNC-BNC Adaptor 5/Pack | \$ 5.25 |
| 40-11-6 | BNC-T Adaptor 5/Pack | \$ 7.35 |
| 40-11-7 | BNC/Double Banana Adaptor | \$ 61.00 |
| 40-11-8 | Stacking Banana Plug - Black 5/Pack | \$ 6.50 |
| 40-11-9 | Stacking Banana Plug - Red 5/Pack | \$ 6.50 |
| 40-12-2 | 12" Stacking Banan/Banana Cord - Black, 5/Pack | \$ 36.75 |
| 40-12-3 | 36" Stacking Banan/Banana Cord - Black, 5/Pack | \$ 42.00 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Devices and Accessories continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|----------------|---|-------------|
| 40-12-4 | 60" Stacking Banan/Banana Cord - Black, 5/Pack | \$ 58.00 |
| 40-12-6 | 12" Stacking Banan/Banana Cord - Red, 5/Pack | \$ 37.00 |
| 40-12-7 | 36" Stacking Banan/Banana Cord - Red, 5/Pack | \$ 42.00 |
| 40-12-8 | 60" Stacking Banan/Banana Cord - Red, 5/Pack | \$ 58.00 |
| 40-13-1 | 12" BNC-BNC Coax (50 Ohms) 5/Pack | \$ 55.00 |
| 40-13-3 | 36" BNC-BNC Coax (50 Ohms) 5/Pack | \$ 60.00 |
| 40-13-9 | BNC Connector only 5/Pack | \$ 15.25 |
| 40-14-1 | Binding Post Bulkhead Conn. - Black, 5/Pack | \$ 9.50 |
| 40-14-2 | Binding Post Bulkhead Conn. - Red, 5/Pack | \$ 9.50 |
| 40-15-05-32-03 | Cable, 4-ft, twisted pair Banana pins to female pins | \$ 68.50 |
| 40-15-1 | Microelectrode Cable 2'/Fem Amph to tin lead (5ea/Pack) | \$ 28.00 |
| 40-40-1 | Xcell 3x4 Probe Only w/2M Cable | \$ 1,345.00 |
| 40-40-3B | XCell 3 Probe (fischer connector) | \$ 795.00 |
| 40-40-7B | Xcell-3 Microelectrode Amp | \$ 2,510.00 |
| 40-40-8B | Xcell 3X4B Amplifier | \$ 6,670.00 |
| 40-41-2 | HiZx8 8-Channel Differential Amplifier | \$ 9,140.00 |
| 40-41-3 | HiZx8 Differential Probe w/2M Cable | \$ 1,440.00 |
| 40-41-5 | Adaptor Cable, DB-25 to 20-pin IDC | \$ 47.00 |
| 40-42-1 | Spike Simulator | \$ 762.00 |
| 40-43-1 | Xcell-3+ Microelectrode Amp | \$ 3,310.00 |
| 40-43-2 | XCell-3+ Probe, 10x | \$ 1,290.00 |
| 40-43-5 | Xcell-3+ Extention Cable | \$ 116.00 |
| 40-46-1 | Spike Signal Enhancer | \$ 1,095.00 |
| 40-50-1B | Regulated Power Supply | \$ 534.00 |
| 40-73-1 | Semi-Chronic Assembly 1 cm | \$ 684.00 |
| 40-73-2 | Semi-Chronic Assembly 2 cm | \$ 810.00 |
| 40-73-2-05 | Semi-Chronic Adaptor, 3/4"Clmn | \$ 451.00 |
| 40-73-7 | Recording Columns, 1 cm dia 5/Pack | \$ 88.00 |
| 40-73-9 | Recording Columns, 2 cm dia | \$ 99.00 |
| 40-73-9-03 | 2.5cm Cap | \$ 39.00 |
| 40-74-0-01 | Column | \$ 138.00 |
| 40-74-0-02 | 1mm Centered Grid | \$ 182.00 |
| 40-74-0-04 | Threaded Cap | \$ 138.00 |
| 40-75-6 | Artifact Suppressor Diff Input | \$ 594.00 |
| 40-77-1B | Microchronic Columns & Caps 3/pack | \$ 242.00 |
| 40-77-2B | Microchronic Recording System | \$ 660.00 |
| 40-77-3 | Recording Barrels 3/Pack | \$ 55.00 |
| 40-77-7 | Electrode Tubing Kit 3/Pack | \$ 159.00 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

p. 9

Devices and Accessories continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|--------------|---|--------------|
| 40-77-8 | Stainless Steel Bone Screws 100/Pack | \$ 0.70 |
| 40-90-2 | 12.5 X 25cm Heating Pad | \$ 136.50 |
| 40-90-2-02 | 25 x 25cm Heating Pad | \$ 156.00 |
| 40-90-2-05 | 10 X 12.5cm Heating Pad | \$ 137.00 |
| 40-90-2-06 | 6.5 X 9.5cm Heating Pad | \$ 137.00 |
| 40-90-2-07 | 5 X 12.5 cm Heating Pad | \$ 132.00 |
| 40-90-3 | Liquid Immersion Heating Rod | \$ 132.00 |
| 40-90-5D-02 | Mini Rectal Thermistor Probe - 1/4" Jack | \$ 158.00 |
| 40-90-6 | Liquid Thermistor Insert | \$ 147.00 |
| 40-90-8D | DC Temperature Controller | \$ 1,642.00 |
| 50-12-10 | Hydraulic Tubing Repair Kit | \$ 126.00 |
| 50-12-1C | Hydraulic Probe Drive | \$ 1,420.00 |
| 50-12-2B | Chronic Adaptor | \$ 116.00 |
| 50-12-3-01 | Trent Wells Conversion | \$ 648.00 |
| 50-12-4 | Acute Mounting Rod w/ Guide Tube Holder | \$ 226.00 |
| 50-12-5 | Replacement Diaphragms 3/Pack | \$ 65.00 |
| 50-12-8-02 | Micrometer Drive | \$ 1,090.00 |
| 50-12-9-02 | Drum Drive | \$ 1,420.00 |
| 50-15-4(IIN) | 4 Electrode Manifold for Multi-Drive | \$ 205.00 |
| 50-15-5(IIN) | 8 Electrode Manifold for Multi-Drive | \$ 290.00 |
| 50-16-1-01 | MCM Controller Module | \$ 2,441.00 |
| 50-16-1-02 | MCU Coupling Unit | \$ 2,940.00 |
| 50-16-1-03 | DRC-01 Digital Remote (Push-Button Control) | \$ 515.00 |
| 50-16-1-04 | ARC-01 Analog Remote (Knob Control) | \$ 610.00 |
| 50-16-1-05 | VRC-01 Virtual Remote (Computer Control) | \$ 179.00 |
| 50-16-1-10 | MCM Accessory Kit | \$ 42.00 |
| 50-16-1-20 | Main Motor Extension Cable 15' | \$ 237.00 |
| 50-16-4-01 | MCM-4 Controller Unit | \$ 6,300.00 |
| 50-16-4-02 | MCU-5 Coupling Unit | \$ 11,200.00 |
| 50-16-4-03 | MCU-4 Coupling Unit | \$ 9,135.00 |
| 50-16-4-10 | MCM-4 Accessory Kit | \$ 103.00 |
| 50-16-4-20 | MCU-5 Accessory Kit | \$ 346.00 |
| 50-16-4-30 | MCU-4 Accessory Kit | \$ 228.00 |
| 50-16-4-40 | 4 Channel Collet Replacement Kit | \$ 132.00 |
| 55-00-1 | Desktop Power Supply | \$ 290.00 |
| 55-00-2 | Power Transfer cable | \$ 32.00 |
| 55-00-3 | Communications cable | \$ 46.50 |
| 55-00-4 | BNC Patch Cable | \$ 26.50 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

p. 10

Devices and Accessories continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|--------------|--|-------------|
| 55-00-5 | Dress Panel For PCF 1" | \$ 20.00 |
| 55-00-6 | Replacement USB Cord | \$ 31.50 |
| 55-00-7 | 2" Blank Panel for SAF | \$ 28.50 |
| 55-00-8 | DB37 M-M Cable 25' | \$ 140.00 |
| 55-00-9 | Male Phono Patch Cable | \$ 31.50 |
| 55-10-0B | APM Acquisition Processor Module | \$ 2,575.00 |
| 55-10-1 | LSG Line Sync Generator | \$ 790.00 |
| 55-11-0 | SAF Rack Frame | \$ 358.00 |
| 55-20-0B | APC-02 Acquisition/Processor Channel | \$ 1,288.00 |
| 55-21-01 | PSC Power Supply Channel | \$ 394.00 |
| 55-21-0B | PCF-02 Powered Channel Frame | \$ 3,255.00 |
| 55-40-0 | PAU 12 Ch Differential Preamp | \$ 1,735.00 |
| 55-40-0-03 | PAU 12 Channel Preamp High Frequency Filtering | \$ 1,735.00 |
| 55-40-1-01 | SPCB-12 12 Channel Switching PCB | \$ 290.00 |
| 55-40-1-04 | SCP Extension Cable 12' | \$ 79.00 |
| 55-40-1-05 | SPS Extension Cable 12' | \$ 79.00 |
| 55-40-2 | SCP Single Channel Preamp | \$ 348.00 |
| 55-40-2-01 | SPS Single Power Supply | \$ 231.00 |
| 55-40-3 | SCM-12 12 Channel Switching Control Module | \$ 1,575.00 |
| 55-40-6(001) | Differential Electrode Leads 2/Pack | \$ 31.50 |
| 55-40-6(006) | Differential Electrode Leads Bare ends | \$ 28.00 |
| 55-50-0B | Experiment Control Module | \$ 8,390.00 |
| 55-50-3 | MSU Monitor Sync Unit see specs | \$ 920.00 |
| 55-50-3-02 | MSU Power Adapter | \$ 63.00 |
| 55-60-0 | PGM Pulse Generator Module Single Channel | \$ 3,950.00 |
| 55-65-0 | ESG Electrical Stimulus Generator | \$ 1,450.00 |
| 55-70-0 | ICM Impedance Conditioning Mod | \$ 1,969.00 |
| 55-70-0-01 | Electrode Conditioning Cable | \$ 142.00 |
| 55-70-0-03 | ICM Alligator Clip Set | \$ 36.75 |
| 55-70-1 | ICM Switching Unit 24 Channel | \$ 499.00 |
| 55-USA | Line Cord, Research 6.5 ft USA For USA use. | \$ 20.00 |
| 60-00-1 | Guide Tube Drive | \$ 1,050.00 |
| 60-00-2 | Guide Tube Drive 45mm Fine Adj | \$ 1,287.00 |
| 60-00-3 | Guide Tube Drive 65mm Fine Adj | \$ 1,383.00 |
| 60-00-4 | GT Drive 45mm w/ Kopf Adaptor | \$ 1,300.00 |
| 60-01-0 | Fixed Rack Base | \$ 184.00 |
| 60-02-0 | Open Chamber XY Drive Base | \$ 1,145.00 |
| 60-03-0 | Closed Chamber XY Drive Base | \$ 1,365.00 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

p. 11

Devices and Accessories continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|------------|---|-------------|
| 60-04-0 | Open Chamber XY Drive | \$ 1,281.00 |
| 60-07-0 | Semi Chronic Adaptor | \$ 184.00 |
| 60-10-0 | Crist Adaptor for XY Drive | \$ 72.00 |
| 61-00-1 | Axial Loading Fixture 23ga tubes | \$ 473.00 |
| 74-20-3B | Diff Amplifier w/Power Supply | \$ 1,777.00 |
| 74-30-1 | Audio Analyzer | \$ 709.00 |
| 74-40-5B | Rate/Interval Monitor | \$ 1,367.00 |
| 74-40-7 | Timed Counter with DC output | \$ 1,633.00 |
| 74-50-8B | Pulsar 6bp Stimulator | \$ 4,226.00 |
| 74-60-1C | S/H Window Discriminator w/Delay Line | \$ 2,750.00 |
| 74-65-6 | bp Optical Isolator with Probe | \$ 2,097.00 |
| 74-65-7 | bp Optical Isolator w/o probe | \$ 1,748.00 |
| 74-77-2 | Glass Microelectrode Carrier | \$ 135.00 |
| 74-77-3 | Metal Microelectrode Carrier | \$ 116.00 |
| 74-77-5 | Metal & Glass Microelectrode Carrier | \$ 176.00 |
| 74-78-1 | Dual Notch Filter | \$ 668.00 |
| 80-10-1 | RT-4 Recording Tank | \$ 1,208.00 |
| 80-10-2 | SC-1 Static Recording Chamber | \$ 210.00 |
| 80-10-3 | Push/Pull Syringe System | \$ 174.00 |
| 80-10-4 | PC-1 Perfusion Chamber, 1 cm | \$ 263.00 |
| 80-10-4-01 | 1cm Sintered Disks | \$ 81.00 |
| 80-10-4-02 | Micro-perfusion Chamber, 1 cm | \$ 263.00 |
| 80-10-6-02 | TC-1 Tissue Chopper Electronic Micrometer | \$ 1,444.00 |
| 80-11-4 | PC-2 Perfusion Chamber, 2 cm | \$ 242.00 |
| 80-11-5 | Replacement Teflon Base for TC-1 | \$ 25.20 |
| 80-11-6 | Replacement Slicer Blades | \$ 2.63 |
| 80-11-7 | PT-3 Pool-type Chamber | \$ 326.00 |
| 80-12-1 | Glass Blade Holder | \$ 135.00 |
| 80-12-2 | Mounting Blocks | \$ 15.00 |
| 80-12-3 | Surface for 2.5" Tray Pedestal | \$ 51.00 |
| 80-12-3-02 | Surface for 3.5" Tray Pedestal | \$ 72.00 |
| 80-12-4 | Replacement Halogen Lamp | \$ 58.00 |
| 80-12-5 | 3.5" (Large) Tray w/Pedestal | \$ 427.00 |
| 80-15-2 | 2.5"(Standard) Tray w/Pedestal | \$ 407.00 |
| 85-00-1 | OTS-5000 Tissue Slicer FHC | \$ 7,718.00 |
| 85-00-2 | Tray for OTS-5000 | \$ 226.00 |
| 85-00-3 | Tray Pedestal | \$ 578.00 |
| 85-00-4 | 2X mag w/ 12" gooseneck FRM18 | \$ 152.00 |

Items on this price list have **not** been cleared for use on humans.

Custom Orders

Please contact us to discuss customizing any of our products to meet your requirements.

International Orders

Taxes, as well as duty, broker, and shipping/handling charges may add significantly to costs on deliveries outside the U.S. For quotes/orders, please contact us to get the name of our agent in your country or work with us directly.

Terms

FOB Bowdoin, ME U.S.A.; Net 30 days, payable in U.S. dollars or by check drawn on U.S. bank, Visa or Mastercard credit cards, wire transfer.

Quantity Discounts

Packaged items -- (PK):
5%, 4-7 pks; 10%, 8-19 pks;
15%, 20 or more pks.

Rush Orders

Standard delivery for FHC products, including microelectrodes, is 2-4 weeks. We will make every effort to ensure that your order arrives in a timely way, and will not charge extra fees for making special trips to after-hours shippers or airports. If you need to receive your order in less than the standard delivery, and the items are not in stock, a 25% rush charge may be added to your order.

Returns

A 20% restocking fee may be charged for returned items. Custom microelectrodes cannot be returned to inventory and, therefore, are not subject to a restocking fee or refund.

Repairs

Repairs are billed at an hourly rate of \$160. A minimum half hour bench charge is applied to all repair work.

Research Products Price List L017W-02 US
Rev. 11/28/11

Devices and Accessories continued

Prices are for each piece; package quantity is listed where appropriate.

| Item | Description/Package Quantity | Price Ea. |
|---------|---|-------------|
| 85-00-5 | Foot Switch | \$ 195.00 |
| 85-00-6 | Sapphire Blade | \$ 431.00 |
| 85-00-7 | Sapphire Blade Clamp for OTS 4000/4500/5000 | \$ 74.00 |
| 85-00-8 | Pedestal Post w/Thumbscrew | \$ 74.00 |
| 85-00-9 | OTS Knurled Blade Clamp Screw | \$ 14.00 |
| 85-12-0 | Tissue Mounting Blocks 1" | \$ 25.00 |
| 85-13-0 | Vise Assembly | \$ 247.00 |
| 85-14-0 | Pedestal Surface | \$ 126.00 |
| 85-20-0 | MRC5000 Medium Refrig. System | \$ 2,284.00 |