


Matrix Bipolar Microelectrode

Parallel configuration using 2 microelectrodes for
Invitro and Invivo Cortical Mapping Single and multi-cell recording and stimulation



Catalog Number	SP (μm)	OD (μm)	R/SF	Material	Ext.	Prot.
30250	50	635	Multi Cell	Pt/Ir	9mm	1mm
30251	115	635	Multi Cell	Pt/Ir	9mm	1mm
30253	200	635	Multi Cell	Pt/Ir	3mm	2mm
30254	250	813	Multi Cell	Pt/Ir	9mm	1mm
30255	250	813	Single Cell	Tungsten	9mm	1mm
30256	500	1245	Single Cell	Tungsten	9mm	1mm