


Monopolar Microelectrode with No Pin Connector

Configured for Use with Guide Tube



Catalog Number	Material	OD/OD Ins. (μm)	Taper	Profile	R/SF	Insulation
30040	Tungsten	200/235	Fine	Thinned	Multi Cell	Epoxy/ite
30041	Tungsten	200/235	Fine	Thinned	Single Cell	Epoxy/ite
30043	Tungsten	125/150	Fine	Thinned	Multi Cell	Epoxy/ite
30044	Tungsten	125/150	Fine	Thinned	Single Cell	Epoxy/ite
30046	Pt/Ir	125/150	Fine	Standard	Multi Cell	Epoxy/ite
30047	Pt/Ir	125/150	Fine	Standard	Single Cell	Epoxy/ite
30049	SS	125/150	Medium	Standard	Multi Cell	Epoxy/ite
30050	SS	125/150	Medium	Standard	Single Cell	Epoxy/ite
30060	Tungsten	100/210	Fine	Standard	Single Cell	Epoxy/PI
30061	Tungsten	100/210	Fine	Standard	User Cond.	Epoxy/PI
30062	Tungsten	75/137	Fine	Standard	Single Cell	Epoxy/PI
30063	Tungsten	75/137	Fine	Standard	User Cond.	Epoxy/PI