


Monopolar Microelectrode with Pin Connector Configured for Acute Recording and Stimulation



Catalog Number	Material	OD/OD Ins. (μm)	Taper	R/SF	Insulation
30020	Tungsten	200/280	Blunt	Multi Cell	EpoxyLite
30021	Tungsten	200/280	Blunt	Single Cell	EpoxyLite
30022	Tungsten	200/280	Blunt	User Cond.	EpoxyLite
30024	Pt/Ir	125/218	Fine	Single Cell	Glass/PI
30025	Pt/Ir	125/218	Fine	Multi Cell	Glass/PI
30026	Pt/Ir	125/218	Medium	Multi Cell	Glass/PI
30030	Tungsten	200/280	Fine	Multi Cell	EpoxyLite
30031	Tungsten	200/280	Fine	Single Cell	EpoxyLite
30032	Tungsten	200/280	Fine	User Cond.	EpoxyLite
30033	SS	250/330	Medium	Single Cell	EpoxyLite
30034	SS	250/330	Medium	Multi Cell	EpoxyLite