



microTargeting™ Solution for Monteris® NeuroBlate® System

Patient-specific stereotactic system that allows a surgeon to quickly and accurately place a Monteris® Mini-Bolt and support the Monteris® Robotic Probe Driver

The microTargeting™ solution for Monteris® NeuroBlate® System is one of the latest technology advances introduced as part of FHC's Connected Care product portfolio and our commitment to working with the medical community in improving patient care. The adapter is designed to act as the interface between the microTargeting Platform and the Monteris NeuroBlate System for MRI-guided neurosurgical ablation procedures.

The adapter is attached to a standard microTargeting Platform via a mounting hub and secured with provided thumbknobs. Once the adapter is correctly oriented in the hub, it enables compatibility with the Monteris Mini-Bolt and Robotic Probe Driver (RPD). Both entry and target mounting hubs attach to the adapter and allow a change in planned trajectory.

microTargeting™ and WayPoint™ are trademarks of FHC, Inc.
Monteris® and NeuroBlate® are registered trademarks of Monteris® Medical

- ✓ **Flexibility:** allows the FHC patient-specific frameless stereotactic platform to be used with the Monteris NeuroBlate System for neurosurgical procedures.
- ✓ **MR Safe:** allows use in MRI environment as well as the OR
- ✓ **Adaptable:** attaches directly to the Monteris NeuroBlate System and to the microTargeting Platform
- ✓ **Steam Sterilizable:** for rapid turn around and ease of use
- ✓ **Convenience:** can be sterilized with the WayPoint™ Sterilization Kit in the SteriSuite Case
- ✓ **Integrated:** for epilepsy cases, reuse scan sets used to implant SEEG depth electrodes



Ordering Information:

MP-KIT-P-NB – Monteris NeuroBlate
microTargeting Platform

Designed for use with:

RPD-01 Monteris Robotic Probe Driver

Additional Items Required for Operation:

C0254 - 3.3mm RPD Adapter

C0255 - 2.2mm RPD Adapter

C0256 - Small Parts Sterilization Basket

C0257 - Ablation Platform Mounting Hardware

L022-03 (Rev C0 2020-03-03)



FHC, Inc.

1201 Main Street
Bowdoin, ME 04287 USA
Fax: +1-207-666-8292
www.fh-co.com



FHC Europe

(TERMOBIT PROD srl)
42A Barbu Vacarescu Str, 3rd Fl
Bucharest 020281 Sector 2
Romania

FHC Latin America

Calle 6 Sur Cra 43 A-200
Edificio LUGO Oficina 1406
Medellín-Colombia